

NOTICE: This document contains references to Varian.
Please note that Varian, Inc. is now part of Agilent
Technologies. For more information, go to
www.agilent.com/chem.



Agilent Technologies

APPLICATION NOTE

Analgesics

Application 2005 LC

Column : Polaris 5 μ C8-Ether, 150 x 4.6mm
Part Nr.:A2030150X046
Eluent : A: H₂O + 0.1% HCOOH
B: MeCN + 0.1% HCOOH
Gradient : 40% to 90%B in 10 min
Flow : 1.0 ml/min
Detection : UV @ 220nm

Peak Identification

1. acetaminophen
2. acetylsalicylic acid
3. salicylic acid
4. ibuprofen

